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Stand der Information: 01/2019 - V6.2

Medizinprodukt: Klasse II b, CE 1250

General information – ActiMaris® wound gel (20g)

1. Product Profile

ActiMaris® wound gel is made of alkaline ionized water with sea salt (3%), NaOCl (0.2%), active oxygen (O₂) as well as lithium-magnesium-sodium-silicate.
ActiMaris® wound gel in interaction with the mucous membrane, skin and wound acts nascent.
Active oxygen, sodium hypochlorite and sea salt have an eliminating effect on typical sore pathogens without damaging the tissue. The decontaminating effect against MRSA / VRE tested positive. The oxidative effect of a singlet oxygen (O₂) and hypertonic saline solution are also responsible for the neutralization of wound odors.
The alkaline pH of ActiMaris® wound gel supported by the provision of OH⁻ ions (as part of a naturally balanced redox system), additionally controls the release of oxygen.
Sea salt contains sodium chloride as its main ingredient (NaCl = 98%) as well as 2% other sea salts.
The increased salinity of ActiMaris® wound gel in turn ensures an optimum hyperosmotic effect so that peri-wound swelling subsides quickly.

2. Fields of application

ActiMaris® wound gel is suitable:
For cleaning, humidification, decontamination, biofilm resolution, swelling and physiological debridement, to create an optimal „micro-environment“ for:

- acute mechanical wounds (skin lacerations, bites, cuts, abrasions, tears and bruises) and postoperative wounds
- chronic wounds (decubitus ulcers, arterial-venous ulcers, diabetic ulcers)
- neurotic, malodorous wounds and ulcerating tumors, even with cavities (Hollow spaces)
- thermal and chemical wounds (burns 1st to 3rd degree)
- entry portals of urological catheters and PEG tubes and drains
- venous ulcers

For the local treatment of skin and mucosa in inflammatory and infectious processes. For the prevention of infection at the mucosa and skin, burns and other wound types. For this, ActiMaris® wound gel is well suited as a wound filler.

3. Instructions for use

ActiMaris® wound gel should be applied at room temperature. Before this, the wound gel can be warmed to body temperature before use. Before application of the wound gel, the wound should be cleaned with ActiMaris® wound irrigation solution (Forte or Sensitive). The wound gel is then liberally applied on or in the wound and covered with an inert compress, gauze etc. The duration and frequency of application should be adapted to the individual wound conditions and should be carried out until all deposits and necrosis are easily removed and the wound is visually clean.

Depending on the wound condition, a bandage change can be done according to indications and phases, every 2 to 3 days. In case of non-infected wounds, necrosis and coatings can be mechanically removed depending on the condition of the wound. In cases of difficult removal of bandages, it is recommended to hydrate the wound dressings with wound irrigation solution ActiMaris® (Forte or Sensitive) for at least 5 minutes until a smooth removal of the bandages is possible without trauma to the wound surface. A primary surgical treatment is recommended with thick, firmly adhering fibrinous or necrotic wound dressings.

For cleaning wounds, ActiMaris® wound gel (see also section 9: Interactions) can be combined with inert dressings or other wound dressings.

4. Tissue tolerance and biocompatibility

ActiMaris® wound gel is suitable for long-term (> 24h) as well as for repeated use during dressing changes, even with irritated skin and mucous membranes. The application is painless and is generally well tolerated by allergy sufferers. Inhibition of granulation and epithelialization is not to be feared. The wound gel is dermatologically assessed as safe.

5. Risks followed by improper use

- Improper use may cause softened wound edges
- Cross-contamination between wounds or patients with incorrect handling of the tube / bottle

6. Side effects

So far, in several years of clinical experience, the following symptoms occasionally occurred

- Slight temporary burning sensation
- If applied to nose area, slight burning and sneezing

7. Contraindications:

Since allergies can never be excluded ActiMaris® irrigation solutions or ActiMaris® wound gel should not be used with known or suspected allergies to any of its ingredients.

8. Restrictions on use

Pregnancy and lactation:

There is no found evidence of genotoxicity, carcinogenicity or embryotoxic effects of the ingredients. Due to the lack of studies and long-term clinical experience with pregnant and breast feeding women, ActiMaris® wound gel should be used in these cases only after an accurate medical assessment.

Infants and newborns:

Due to insufficient experience on long-term exposure, ActiMaris® wound gel should only be applied after careful medical benefit-risk assessment.

9. Interactions

ActiMaris® wound gel acts as a bio-physical compound due to the content of ROS substances (O₂ + OCl⁻), energy and oxidizing agents, which is why a simultaneous combination with other local, active substances and / or interactive wound dressings is in principle unnecessary. Contact with non-precious metals should be avoided. The wound can be covered with modern wound dressings after applying ActiMaris® wound gel.

On the other hand, a good compatibility with modern wound dressings such as silver dressings has been shown in clinical trials. Likewise, there is ample clinical evidence of good tolerability when combined with inert, non-active bandages and wound dressings such as compresses, gauze, pads, sponges, gels, hydro fibers, alginates, hydrocolloids etc. Concomitant use of ActiMaris® wound irrigation solution with the ActiMaris® wound gel is possible and useful.

10. General safety instructions

ActiMaris® wound gel is suitable only for external use of wounds, skin or mucosa. Do not use for infusion or injection! Do not swallow! Use only closed and undamaged containers. Keep out of reach of children! It is possible for colored fabrics to fade after contact with ActiMaris® products.

11. Composition and technical information

Sea salt, water, lithium-magnesium-sodium-silicate, sodium oxchlorite.
Clear, colorless, alkaline and hypertonic gel with characteristic smell; contains 3% sea salt (Sal Maris) and 0,2% NaOCl.

12. Storage specification, duration, expiry date

ActiMaris® wound gel has to be stored at room temperature between +5° and +25°C. Protect from direct sunlight. Shelf life is 30 months. Expiry date is on the packaging. After opening the tube, the expiration period is 3 months; always close tube tightly after use.

ActiMaris AG Manufacturer:

ActiMaris AG
Sandgrube 29
9050 Appenzell
info@actimaris.com
Tel.: +41 71 505 75 25

- LOT Lot number
- REF Order number
- Follow instructions for use
- Use by: year and month
- Storage (temperature)
- Protect from direct sunlight
- Manufacturer

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ActiMaris AG **Fababricant:**